Sterilization Container Systems



safe, secure, sterilized

Protecting Your Instruments Every Time

5 Different Models





STERILIZATION CONTAINER SYSTEMS

Hermann Müller window model sterilization container system is one of its patented products. What makes this model different from the others is the window on the lid which its name comes from. This window of made of high performance PPSU transparent material which is resistant against high temperature, impacts, scratches, deformation, discoloration and chemicals and it opens users a way to visually check inside of the sterilization container without opening the lid. The window model lids assures great tightness due to its special assembly method and the gasket system in it. As there is no filtering mechanism on the lid, window model sterilization container boxes strictly must be bottom perforated. Either single use paper or 2000 cycles PTFE Teflon Filters can be used in the box. In order to allow the air flow and reach

to perfect steam penetration the filtering mechanism of this model must have enough space. For that purpose, this sterilization container model is supplied with a wire base inside the container as standard. Except for this container model, all Hermann Müller boxes and lids are compatible with each other.

Hermann Müller window model sterilization container systems are manufactured as five different length, width and height. All models and sizes have specific baskets and accessories which can be viewed in following parts. Modular auto lock mechanism (thermo lock) promising sterility assurance and indicator window providing visibility of the indicator card result by visual control without opening the lid are also available for this sterilization container model.

Box Types



Box (With Wire Base)



Standard Heights











111 mm

131 mm

153 mm

205 mm

260 mm





Dimensions







1/1 Full Size

Page 60 -61

3/4 Quarter Size

Page 62 -63

1/2 Half Size

Page 64 -65

STERILIZATION CONTAINER SYSTEMS

1/1 Full Size

Box Types



Box (With Wire Base)







Box (With Wire Base)



Code	Outside Dimension	Inside Dimension	Weight
K6-100D	600 x 272 x 111 mm	543 x 258 x 66 mm	2935 g
K6-101D	600 x 272 x 131 mm	543 x 258 x 82 mm	3095 g
K6-102D	600 x 272 x 153 mm	543 x 258 x 104 mm	3255 g
K6-103D	600 x 272 x 197 mm	543 x 258 x 156 mm	3710 g
K6-104D	600 x 272 x 250 mm	543 x 258 x 201 mm	4090 g

Window Lids



Color	Code	Color	Code
Grey	K6-120D	Red	K6-124D
Yellow	K6-121D	■■ Black	K6-125D
Green	K6-122D	Purple	K6-126D
Blue	K6-123D	Bronze	K6-127D
Dimension 6	602 x 285	mm Weight	1305 g

STERILIZATION CONTAINER SYSTEMS

3/4 Quarter Size

Box Types



Box (With Wire Base)







Box (With Wire Base)



Code	Outside Dimension	Inside Dimension	Weight
KC 300D	475 v 272 v 444 mm	422 v 258 v 66 mm	2260 ~
K6-300D	475 x 272 x 111 mm	422 X 258 X 66 IIIIII	2360 g
K6-301D	475 x 272 x 131 mm	422 x 258 x 82 mm	2450 g
K6-302D	475 x 272 x 153 mm	422 x 258 x 104 mm	2585 g
K6-303D	475 x 272 x 200 mm	422 x 258 x 151 mm	3040 g

Window Lids



Color	Code	Color	Code
Grey	K6-340D	Red	K6-344D
Yellow	K6-341D	■■ Black	K6-345D
Green	K6-342D	Purple	K6-346D
Blue	K6-343D	Bronze	K6-347D

Dimension 477 x 285 mm **Weight** 1020 g

STERILIZATION CONTAINER SYSTEMS

1/2 Half Size

Box Types



Box (With Wire Base)







Box (With Wire Base)



Code	Outside Dimension	Inside Dimension	Weight
K6-600D	310 x 272 x 111 mm	253 x 258 x 66 mm	1765 g
K6-601D	310 x 272 x 131 mm	253 x 258 x 82 mm	1825 g
K6-602D	310 x 272 x 153 mm	253 x 258 x 104 mm	1935 g
K6-603D	310 x 272 x 197 mm	253 x 258 x 156 mm	2185 g
K6-604D	310 x 272 x 250 mm	253 x 258 x 201 mm	2435 g

Window Lids



Color	Code	Color	Code
Grey	K6-630D	Red	K6-634D
Yellow	K6-631D	Black	K6-635D
Green	K6-632D	Purple	K6-636D
Blue	K6-633D	Bronze	K6-637D
Dimension 3	312 x 280 n	nm Weight	736 g



Headquarter

Hermann Müller Medizintechnik GmbH Dr. Karl-Storz-Strasse 3 78532 • Tuttlingen • Germany

T:+49 7461 4125 F:+49 7461 79263

Factory

Kerimbeyosb Mah. Yaşardoğu Cad. No: 76/1 55330, Tekkeköy/Samsun/Türkiye

www.hermannmuller.com info@hermannmuller.com

All rights reserved.No part of this catalog may be reproduced in any from without the written permission from Hermann Müller